

ROWHEELS® **REV** Order Form

PROVIDER INFORMATION		RTS INFORMATION	
Provider:		Name:	
Provider POC:		Email:	
Account #:		Phone #:	
Email:		Address:	
Phone #:		City:	State: ZIP:
Fax #:		Ship to information	
Address:		Company (if applicable):	
City:	State: ZIP:	Name:	
Client:		Phone #:	
Rowheels information		Address:	
Ship date:		City:	State: ZIP:
Order #:		Ship via:	
User Information		Wheelchair Information	
Weight:	Height:	Wheelchair Make:	Model:
Med. Diagnosis:		Camber Angle:	Other:

Rowheels are covered through HCPSC billing code E2227

Please email or fax to Rowheels customer service at orders@rowheels.com or (608) 719-0010



Key 1: MODEL

MODEL	KEY CODE	KEY INPUT	
REV-LX (Low gear).....	RWLX		\$5,400 (E2227)
REV-HX(High gear).....	RWHX		\$5,400 (E2227)



The Right Way to Roll.™

www.rowheels.com

8001 Terrace Ave Suite #204
Middleton, WI 53562

PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com

ROWHEELS® **REV** Order Form

Key 2: COLOR SCHEME

COLOR SCHEME	KEY CODE	KEY INPUT
Red/Black.....	0.....STD	
Black/Black.....	1.....STD	
White/Black.....	2.....STD	
Red, Black or White/Custom	C.....RFQ (Custom hydrographic on tri-adapter)	



Red/Black



Black/Black



White/Black

Key 3: TORQUE ARM

TORQUE ARM	KEY CODE	KEY INPUT
8" Torque Arm.....	8.....STD	
10" Torque Arm.....	10.....STD	

Key 4: HANDRIM

HANDRIM	KEY CODE	KEY INPUT
Out-Front Natural-Fit.....	N.....\$335	
Out-Front Natural-Fit LT.....	NL.....\$365	
Out-Front SURGE.....	S.....\$425	
Out-Front SURGE LT.....	SL.....\$475	
Q-Grip.....	Q.....\$315 (not available on 22" or 26" wheels)	
RibGrips.....	G.....\$260 (only available on 24" wheels & in blue, pink, black and red)	
Projection-8 (45°).....	P.....\$375	



The **Natural-Fit** has an oval profile that increases gripping surface area by 80%. The optional thumb grip component is not available for use on Rowheels.

The **Surge** has the same oval profile along with a rubber strip along the top surface of the oval for increased friction while leaving the side surface smooth for braking.

Q-Grips are ideal for anyone with limited hand function or anyone that prefers the feel & functionality of a fully-coated handrim. Research shows the Q-Grip to be 59.1% more grippy than a vinyl handrim.

RibGrips are ideal for those looking to push their wheelchair with more ease and comfort. The soft rib discs benefit users with poor hand/finger grip.



The Right Way to Roll.™

www.rowheels.com

8001 Terrace Ave Suite #204
Middleton, WI 53562

PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com

ROWHEELS® **REV** Order Form

Key 5: AXLE LENGTH

AXLE LENGTH	KEY CODE	KEY INPUT
Length 0 (2.0").....	0.....STD	<i>*Length 0 fits most Tilite, Kimobility, Pride frames</i>
Length 1 (2.1").....	1.....STD	<i>*Length 1 fits most Colours, Motion Composites frames</i>
Length 2 (2.2").....	2.....STD	<i>*Length 2 fits most Quickie frames</i>
Length 3 (2.3").....	3.....STD	<i>*Length 3 fits Top End Crossfire/T6 frames</i>
Length 4 (2.4").....	4.....STD	<i>*Length 4 fits Invacare ProSpin/Patriot Plus, Top End All-Terrain</i>
Length 5 (2.5").....	5.....STD	

Key 6: TIRE

TIRES	KEY CODE	KEY INPUT
No tire.....	N.....NA	
Schwalbe RightRun (Black/Grey)...R.....	STD	
Schwalbe Marathon Plus (Black)...M.....	\$130	
Shox Solid (Black).....S.....	\$59	
Kenda Nevegal.....X.....	\$78 (not available on 22" & 26" wheels)	



The Schwalbe **RightRun** tires incorporate their most effective SmartGuard puncture protection layer, hand-friendly 2Grip side walls, and attractive two color compound. Robust Black'n Roll tread center compound that is non-marking



The Schwalbe **Marathon Plus** is the benchmark for puncture-proof pneumatic tires on the wheelchair. The all-new version comes with improved rolling resistance along with hand-friendly 2Grip side walls, and the robust Black'n Roll tread center compound.



SHOX Solid Wheelchair Tires are airless so you don't have to worry about punctures and are a firm 110psi equivalent.



The Kenda **Nevegal** is the optimal tire for an assortment of ever-changing terrain. It is highly adaptable and excels under various conditions, from loose rock, tall grass, muddy roots, single track and high speed descents, and many more.

Key 7: WHEEL SIZE

WHEEL SIZE	KEY CODE	KEY INPUT
22" (489mm).....	2.....STD	
24" (540mm).....	4.....STD	
25" (559mm).....	5.....STD	
26" (590mm).....	6.....STD	



The Right Way to Roll.™

www.rowheels.com

8001 Terrace Ave Suite #204
Middleton, WI 53562

PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com

ROWHEELS® **REV** Order Form

KEY 8: SPOKE GUARD

SPOKE GUARD	KEY CODE	KEY INPUT
Carbon Fiber*	S (STD)	
Clear	C (STD)	
Custom	Z (RFQ)	

*Guard is covered in a vinyl carbon fiber wrap, Not actual carbon fiber



FINAL ORDER KEY CODE

Enter, in order, the seven option keys selected on the previous page here:

KEY#1 KEY#2 KEY#3 KEY#4 KEY#5 KEY#6 KEY#7 KEY#8

Sample Final Order Key code - **RWLX-0-8-N-1-R-4**


Description - REV 1, Red/Black, 8" Torque Arm, Natural-Fit, Length 1, Right Run Tire, 24"

REV¹ WHEELCHAIR ADAPTER INFO

In order for REV wheels to function properly, our universal frame clamp latches must be installed on the wheelchair's frame. If the chair to be used has a tubular frame with a diameter between **0.75" (19.05mm) and 1.375" (34.93mm) or is a Quickie Q7**, select "Compatible Frame", otherwise select "Not Compatible" and provide us with the wheelchair make, model and frame dimensions in the "Notes" section below. A decision regarding whether a custom adapter can be designed will be made. An RFQ can then be submitted for the cost of the custom adapter.

Compatible Frame

Not Compatible



NOTES:

The Right Way to Roll.™

www.rowheels.com

8001 Terrace Ave Suite #204

Middleton, WI 53562

PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com

